



Natural gas analysis

BIOMETHANE
REFINERY GASES
NATURAL GASES
INDUSTRIAL EMISSIONS AND VOCS
INDUSTRIAL AND SPECIALITY GASES

APPLICATIONS

Certified measurement of Gross Calorific Value of Natural gas processes.

Through metrological certification (OIML R140, Welmec 7.2 et ISO 6976), our device is able to provide the Calorific Value that leads to natural gas tariffication.

2 modules are mandatory for this 2 min analysis.


PixLPro software gives you direct overview of the measurement.

On the ChromPix® or ChromEx 400 device, you can add two optionnal modules to analyze sulfur compounds (H₂S and odorizant - TBM or THT)

SYSTEMS



ChromPix2®



TwinPix®



ChromEx200/400®



SAMPLE

Typical natural gas sample

CONCLUSIONS

APIX systems are pending certified for the measurement of Gross Heating Value of Natural Gas.

Two analytical modules (PPU10 and PDMS10) are mandatory for performing this certified measurement.

System configuration can be completed with two other modules for providing analysis of supplementary compounds (sulfur compounds...)


MODULE A

Reference	MK10-TCD-2μL-PPU10-PPU1-F1
Detector	TCD
Column	PPU (Pora-Plot U) 10m
Internal Diameter	0.25mm
Phase Thickness	12μm
Precolumn	PPU (Pora-Plot U) 1m
Backflush	✓
Sample Loop	2μL
Regeneration	✗

METHOD

Carrier Gas	Helium
Carrier Gas Pressure (max)	36.2 psi - 2.5 bar
Detector T°C	70°C
Column T°C	70°C
Column Pressure	0.8 bar
Sample Loop T°C	70°C
Sample Loop Pressure	0.5 bar
Injection Time	10s
Analysis Time	120s

CHROMATOGRAM



RESULTS

	LOD	LOQ	RSD
1 N ₂ +O ₂	2ppm	6ppm	0.9% (0.6%)
2 CH ₄	1%	3%	0.05% (82.81%)
3 CO ₂	2ppm	6ppm	0.25% (0.29%)
4 C ₂ H ₆	2ppm	6ppm	0.25% (11.81%)


MODULE B

Reference	MK6-TCD-2μL-PDMS10-F1
Detector	TCD
Column	PDMS (100% Methyl Polysiloxane) 10m
Internal Diameter	0.15mm
Phase Thickness	1.2μm
Precolumn	None
Backflush	✗
Sample Loop	2μL
Regeneration	✗

METHOD

Carrier Gas	Helium
Carrier Gas Pressure (max)	36.2 psi - 2.5 bar
Detector T°C	70°C
Column T°C	60°C
Column Pressure	1.2 bar
Sample Loop T°C	60°C
Sample Loop Pressure	0.5 bar
Injection Time	1s
Analysis Time	100s

CHROMATOGRAM



RESULTS

	LOD	LOQ	RSD
1 Injection			
2 C3	1ppm	3ppm	0.1% (7%)
3 iC4	1ppm	3ppm	0.2% (0.1%)
4 nC4	1ppm	3ppm	0.1% (2%)
5 neoC5	1ppm	3ppm	0.3% (0.35%)
6 iC5	1ppm	3ppm	0.2% (0.1%)
7 nC5	1ppm	3ppm	0.4% (0.15%)
8 nC6	1ppm	3ppm	0.6% (500ppm)


MODULE C

Reference	MK10-TCD-20µL-PPU10-PPU1-F1
Detector	TCD
Column	PPU (Pora-Plot U) 10m
Internal Diameter	0.25mm
Phase Thickness	12µm
Precolumn	PPU (Pora-Plot U) 1m
Backflush	✓
Sample Loop	20µL
Regeneration	✗

METHOD

Carrier Gas	Helium
Carrier Gas Pressure (max)	36.2 psi - 2.5 bar
Detector T°C	80°C
Column T°C	105°C
Column Pressure	1.5 bar
Sample Loop T°C	105°C
Sample Loop Pressure	0.5 bar
Injection Time	13s
Analysis Time	80s

CHROMATOGRAM



RESULTS

		LOD	LOQ	RSD
1	Injection			
2	H ₂ S	1ppm	3ppm	10% (3ppm)
3	COS	1ppm	3ppm	10% (3ppm)


MODULE D

Reference	MK6-TCD-20µL-PDMSP12-F2
Detector	TCD
Column	PDMSP (20 % Diphenyl - 80 % Methylpolysiloxane) 10m
Internal Diameter	0.25mm
Phase Thickness	1µm
Precolumn	None
Backflush	✗
Sample Loop	20µL
Regeneration	✗

METHOD

Carrier Gas	Helium
Carrier Gas Pressure (max)	36.2 psi - 2.5 bar
Detector T°C	60°C
Column T°C	60°C
Column Pressure	1.5 bar
Sample Loop T°C	60°C
Sample Loop Pressure	0.5 bar
Injection Time	1s
Analysis Time	200s

CHROMATOGRAM



RESULTS

		LOD	LOQ	RSD
1	Injection			
2	TBM	1ppm	2ppm	<10% (2ppm)
3	THT	1ppm	2ppm	<10% (2ppm)